


## Preparation of Competent Cells

### Competent cells:

*E. coli* are more likely to incorporate foreign DNA if their cell walls are altered so that DNA can pass through more easily. Such cells are said to be “**competent**”. Cells are made competent by a process that uses **calcium chloride (CaCl<sub>2</sub>)** and **heat shock**. Cells that are undergoing very rapid growth (@ Log phase) are made competent more easily than cells in other stages of growth.

[http://www.phschool.com/science/biology\\_place/labbench/lab6/competen.html](http://www.phschool.com/science/biology_place/labbench/lab6/competen.html)



### Procedure:

1. Streak from DB 3.0 stock onto LB plate. Incubate overnight at 37°C, 200 rpm
2. Pick a single colony and inoculate a 2 ml culture of SOB. Incubate overnight at 37°C
3. Add 20 ul overnight culture to **50 ml** SOB. Grow until OD<sub>600</sub>=0.6-0.8 (37°C, 200 rpm)  
*Protocol from Yang Lab: add 2 ml overnight culture to 500 ml SOB, and shake at 23°C; 80 ml TfbI resuspension; 30-50 ml TfbII resuspension*
4. Put the flasks on ice for **10 min** to arrest the *E. coli* growth. (*From now on, keep everything as cold as possible and work on ice*)
5. Divide the bacterial solution into 16 sterile 1.5 ml Eppendorf tube; add **1 ml** bacterial solution into each tube; centrifuge at **6,000 rpm for 5 min at 4°C**; repeat 3 times to pellet all 50 ml SOB  
*Protocol from Yang Lab: centrifuge at 4,000 rpm for 10 min at 4°C*
6. Discard supernatant carefully, and for each tube, *gently* resuspend cell pellet in **500 ul TfbI** while keep on ice
7. Incubate on ice for **15 min** (*Protocol from Yang Lab: 30 min*)
8. Pellet cells at **4°C for 10 min at 4000 rpm**
9. For each tube, resuspend cell pellet in 312 ul TfbII while keep on ice
10. Incubate on ice for **15 min** (*Protocol from Yang Lab: 30 min*)
11. Aliquot **100 ul** into new sterile 1.5 ml Eppendorf tube
12. Flash freeze in liquid nitrogen and store at -80°C

**SOB solution**

	For 1 L	For 500 ml	Final [c]
Tryptone	20.0 g	10.0 g	
Yeast extract	5.0 g	2.5 g	
NaCl (Sodium chloride)	0.6 g	0.3 g	10 mM
KCl (Potassium chloride)	0.5 g	0.25 g	6.7 mM
MgSO <sub>4</sub> (Magnesium sulfate)	1.2 g	0.6 g	10 mM
MgCl <sub>2</sub> •6H <sub>2</sub> O (Magnesium chloride hexahydrate)	2.0 g	1.0 g	10 mM
pH 7.0, NaOH			

\* sterilize the solution through autoclaving

\* Tryptone (胰蛋白胨) is the assortment of peptides formed by the digestion of casein by the protease trypsin. It provides a source of amino acids for the growing bacteria.

\* Yeast extract (酵母粉) is a complex and widely used hydrolysate of yeasts. It provides nitrogenous compounds, carbon, sulfur, trace nutrients, vitamin B complex and other important growth factors, which are essential for the growth of diverse microorganisms.

**Tfbl (transformation buffer I)**

	For 400 ml	For 200 ml	Final [c]
KAc (Potassium acetate)	1.1778 g	0.5889 g	30 mM
MnCl <sub>2</sub> (Manganese(II) chloride)	3.958 g	1.979 g	79 mM
KCl (Potassium chloride)	2.984 g	1.492 g	100 mM
CaCl <sub>2</sub> (Calcium chloride)	0.588 g	0.294 g	13 mM
Glycerol	60 ml	30 ml	15% (v/v)
pH 5.8, Acetic acid			


\* sterilize the solution using syringe filter with a 0.22 µm pore size

**Tfbll (transformation buffer II)**

	For 200 ml	For 100 ml	Final [c]
MOPS	0.418 g	0.209 g	10 mM
CaCl <sub>2</sub> (Calcium chloride)	2.205 g	1.1025 g	100 mM
KCl (Potassium chloride)	0.1492 g	0.0746 g	10 mM
Glycerol	30 ml	15 ml	15% (v/v)
pH 6.5, NaOH			

\* sterilize the solution using syringe filter with a 0.22 µm pore size

\* MOPS is frequently used as a buffering agent in biology and biochemistry.



So B	1 liter	500 mL	250 mL
Tryptone	20.0 g	10 g	5
Yeast Extract	5.0 g	2.5 g	1.25
NaCl (sodium chloride)	0.6 g	0.3 g	0.15
KCl (potassium chloride)	0.5 g	0.25 g	0.125
MgSO <sub>4</sub> (10mm) (Magnesium Sulfate Anhydrous)	1.2 g	0.6 g	0.3
MgCl <sub>2</sub> · 6H <sub>2</sub> O (10mm) (Magnesium chloride)	2.0 g	1.0 g	0.5
pH 7.0 ↓ for MgCl <sub>2</sub> 0.94 g/L			
Chemical method buffer			
TfBI	400 mL	200 mL	
KAC	1.178	0.5889 g	potassium acetate
MnCl <sub>2</sub>	3.958	1.979 g	Manganese dichloride
KCl	2.984	1.492 g	potassium chloride
CaCl <sub>2</sub>	0.588	0.294 g	Calcium chloride
glycerol	60 mL	30 mL	
pH 5-8 acetic acid			
TfBI	200 mL	100 mL	200 mL
Na-Mops	0.418	0.209 g	mops $\frac{186}{219.26} \times 0.418 = 0.3715$
CaCl <sub>2</sub>	2.205	1.1025 g	
KCl	0.1492	0.0746 g	
glycerol	30 mL	15 mL	
pH 6-5 NaOH			
过滤除菌. 0.22 μm filter.			

Start overnight LB culture (2-5 mL)  
 Make sure have 2-50mL LB flasks ready for tomorrow

Electrical Competent Cell

50B. (1L)	(1L)	500 mL
tryptone	20g	10g
Yeast extract	5g	2.5g
NaCl	0.584g	0.292g
KCl	0.1864g	0.0932g

10% Glycerol. ~~125~~ 125g/mL.

1 L add 100mL glycerol = 125g.