


## ***Agrobacterium*-mediated *Tragopogon* Transformation**



### **Protocol**

#### ***Agrobacterium* Culture**

- Streak the *Agrobacterium tumefaciens* onto a LB + Kan + Rif plate. Grow at 28 °C for 2 days.
  - pCambia 1300 vector has Kan- and Hyg-resistance genes; *Agrobacterium* strain EHA105 has Rif-resistance gene
  - To knock out the *PDS* gene, bacteria stock 093 (1300-AtCas9-GmUbi-GFP-AtU6-1-gTraPDS2-AtU6-29-gTraPDS1) will be used
- *Agrobacterium* colonies are isolated and inoculated into 20 ml liquid Ty medium with AS and appropriate antibiotics (in an autoclaved flask). Shake (dark condition) at 100 rpm, 28 °C for ~24 h.
  - Avoid large lumps of bacteria; brush the bacteria against the inside surface of the flask
  - The streaked *Agrobacterium* plate can be used for 1-2 weeks, but no longer than that (e.g., 1 month)

#### **Infection Medium Preparation**

- Resuspend *Agrobacterium* pellet in liquid AAM infection medium (containing AS).
  - Add 1 ml *Agrobacterium* culture to a 1.5 ml microcentrifuge tube
  - Centrifuge at 2,000 rpm for 10 min
  - Remove the Ty medium, and re-suspend the *Agrobacterium* pellet with 1 ml AAM medium

Updated on 11/10/2021

- Adjust the *Agrobacterium* cell density to  $OD_{600} = 0.1$  (blank: AAM + AS without *Agrobacterium*)

## Explant Preparation

- Place wounded leaf segments (i.e., explants) on **M5 medium** in dark for 1-2 days at room temperature.
  - Wounded-side down; shallow wounds act as sites for *Agrobacterium* infection

## Explant Infection

- Immerse explants in a Petri dish containing *Agrobacterium* suspension culture for 30 min at room temperature.
  - Transfer ~20 ml AAM infection medium (containing *Agrobacterium*) into 100 ml X 15 ml Petri dish
  - Swirl the *Agrobacterium* suspension culture (with the explants) every 10 min
  - Blot explants completely on a sterile filter paper; this step usually takes ~20 min
- Transfer infected explants to **Cocultivation medium** covered with autoclaved filter paper (preventing rampant growth of *Agrobacterium*, and therefore, avoid tissue necrosis/browning) and put the plates in dark for 2-3 days at room temperature.

*Callus induction, shooting and rooting are carried out at 25 °C with 14 h photoperiod in an incubator*

## Callus induction

- Transfer the explants to **Selective M5 medium** for callus induction; it usually takes 3 weeks.

## Shooting

- Transfer the callus to **Selective COCO-1 medium** for 2 weeks for shoot induction.
- Transfer the callus to **Selective WPM medium** for shoot formation; shoots will start forming after 2 weeks.

Updated on 11/10/2021

## Rooting

- Excise regenerated shoots from the callus and transfer the shoots to [COCO-R medium](#).
- Transfer the rooted shoots to soil.
  - Transfer rooted shoots to soil inside a tray with a lid on. Put the tray in the growth chamber for 2 weeks. Repot the plants and put the pots inside the greenhouse

## Timeline

Day	<i>Agrobacterium</i> preparation	Explant preparation
1	Streak the stock on plate in the morning; take ~2 days to grow	
2		Place leaf segments on M5 medium in dark for 1-2 days at RT
3	Pick colonies and shake in Ty medium for 24 hours	
4	Prepare the AAM with <i>Agrobacterium</i> and transform the explants; cocultivation of explants and <i>Agrobacterium</i> in dark for 2-3 days at RT	
5		
6	Transfer explants to selective M5 medium	
27	Transfer callus to selective COCO-1 medium	
41	Transfer callus to selective WPM medium for shooting	
62	Transfer shoots to COCO-R medium for rooting	
83	Move regenerated seedlings to soil	

## Abbreviations

Abbreviation	Description
AS	Acetosyringone (乙酰丁香酮): a signal (phenolic compound) exuded by wounded eudicot plant cells and involved in plant-pathogen recognition; AS is recognized by a receptor encoded by the <i>virA</i> gene, which facilitates the processing and transfer of T-DNA from <i>Agrobacterium</i> to the nuclear genome of the plant host; AS, a <i>vir</i> gene inducer, allows higher transformation efficiency
Cefo	Cefotaxime (头孢噻肟): often used in transformation of plant species for the elimination of <i>Agrobacterium tumefaciens</i>
Hyg	Hygromycin B (潮霉素 B): an antibiotic produced by <i>Streptomyces hygroscopicus</i> ; active against both prokaryotic and eukaryotic cells by inhibiting polypeptide synthesis
Kan	Kanamycin (卡纳霉素): works by interfering with protein synthesis – it binds to the 30S subunit of the bacterial ribosome
LB	Lysogeny broth
Rif	Rifampicin (利福平): an antibiotic works by decreasing the production of RNA by bacteria
Tim	Timentin (特美汀): used most commonly in the regeneration medium for elimination of the <i>Agrobacterium</i> post-transformation

## Medium

### Ty medium (liquid)

	For 1 L	For 200 ml
Tryptone	5 g	1 g
Yeast extract	3 g	0.6 g
pH 5.5 (use 1 M HCl)		
Before use add:		
<ul style="list-style-type: none"> <li>• Rif (1,000 X; 50 mg/ml); final [c] = 50 mg/L</li> <li>• Kan (1,000 X; 50 mg/ml); final [c] = 50 mg/L</li> <li>• AS (1,000 X; 100 mM); final [c] = 0.1 mM</li> </ul>		

### AAM infection medium (liquid)

	For 1 L	For 500 ml
AA macro (10 X)	100 ml	50 ml
AA micro (1,000 X)	1 ml	0.5 ml
AA amino acid (100 X)	10 ml	5 ml
MS vitamins (100 X)	10 ml	5 ml
Fe <sub>2</sub> EDTA (100 X)	10 ml	5 ml
Casamino acids	0.5 g	0.25 g
Sucrose	68.5 g	34.25 g
Glucose (0.5 g/ml)	72 ml	36 ml
pH:		
<ul style="list-style-type: none"> <li>• 5.33 before autoclaving (using KOH); add filter-sterilized glucose; the final pH 5.15 (should be around 5.2)</li> </ul>		
Before use add:		
<ul style="list-style-type: none"> <li>• AS (1,000 X; 100 mM); final [c] = 0.1 mM</li> </ul>		

### M5 medium

	For 1 L	For 400 ml
MS salt	4.44 g	1.776 g
Sucrose	30 g	12 g
BAP (1 mg/ml)	1 mg	0.4 mg
NAA (1 mg/ml)	0.5 mg	0.2 mg
Gelrite	5 g	2 g
pH: set to 5.8		

### Cocultivation medium

	For 1 L	For 400 ml
MS salt	4.44 g	1.776 g
Sucrose	30 g	12 g
BAP (1 mg/ml)	1 mg	0.4 mg
NAA (1 mg/ml)	0.5 mg	0.2 mg
Gelrite	5 g	2 g
pH: set to 5.35; will reach 5.2 (the desired pH) after autoclaving Add AS (1,000 X; 100 mM) after autoclaving; final [c] = 0.1 mM		

### Selective M5 medium

	For 1 L	For 400 ml
MS salt	4.44 g	1.776 g
Sucrose	30 g	12 g
BAP (1 mg/ml)	1 mg	0.4 mg
NAA (1 mg/ml)	0.5 mg	0.2 mg
Gelrite	5 g	2 g
pH: set to 5.8 After autoclaving, add: <ul style="list-style-type: none"> <li>• Tim (500 X; 200 mg/ml); final [c] = 300 mg/L</li> <li>• Hyg (50 mg/ml stock) <ul style="list-style-type: none"> <li>○ 15 mg/L Hyg for <i>T. porrifolius</i>; add 300 µl to 1 L medium</li> </ul> </li> </ul>		

### Selective COCO-1 medium

	For 1 L	For 400 ml
Woody Plant Medium	2.41 g	0.964 g
Sucrose	30 g	12 g
BAP (1 mg/ml)	2.5 mg	1 mg
Agargellan	8.5 g	3.4 g
Volume	950 ml (+ 50 ml coco)	380 ml (+ 20 ml coco)
pH: set to 5.7 After autoclaving, add: <ul style="list-style-type: none"> <li>• filter-sterilized coconut water (final 5%)</li> <li>• Tim (500 X; 200 mg/ml); final [c] = 300 mg/L</li> <li>• Hyg (50 mg/ml stock) <ul style="list-style-type: none"> <li>○ 15 mg/L Hyg for <i>T. porrifolius</i>; add 300 µl to 1 L medium</li> </ul> </li> </ul>		

Updated on 11/10/2021

### Selective WPM medium

	For 1 L	For 400 ml
Woody Plant Medium	2.41 g	0.964 g
Sucrose	30 g	12 g
Agargellan	5 g	2 g
pH: set to 5.7 After autoclaving, add: <ul style="list-style-type: none"><li>• Tim (500 X; 200 mg/ml); final [c] = 300 mg/L</li><li>• Hyg (50 mg/ml stock)<ul style="list-style-type: none"><li>○ 15 mg/L Hyg for <i>T. porrifolius</i>; add 300 µl to 1 L medium</li></ul></li></ul>		

### COCO-R medium

	For 1 L	For 400 ml
Woody Plant Medium	2.41 g	0.964 g
Sucrose	30 g	12 g
IBA (1 mg/ml)	1.5 mg	0.6 mg
Agargellan	8.5 g	3.4 g
pH: set to 5.7		

### AA Macro 10X solution

	For 1 L	For 500 ml
KCl	29.5 g	14.75 g
MgSO <sub>4</sub> ·7H <sub>2</sub> O	2.5 g	1.25 g
CaCl <sub>2</sub> ·2H <sub>2</sub> O (or CaCl <sub>2</sub> )	1.5 g (or 1.13 g)	0.75 g (or 0.565)
NaH <sub>2</sub> PO <sub>4</sub> ·H <sub>2</sub> O	1.5 g	0.75 g
Store at 4 °C		

### AA Micro 1000X solution

	For 1 L	For 500 ml
MnSO <sub>4</sub> ·4H <sub>2</sub> O	10 g	5 g
H <sub>3</sub> BO <sub>3</sub>	3 g	1.5 g
ZnSO <sub>4</sub> ·7H <sub>2</sub> O	2 g	1 g
KI	750 mg	375 mg
Na <sub>2</sub> MoO <sub>4</sub> ·2H <sub>2</sub> O	250 mg	125 mg
CuSO <sub>4</sub> ·5H <sub>2</sub> O (or CuSO <sub>4</sub> )	25 mg (or 15.98 mg)	12.5 mg (or 7.99 mg)
CoCl <sub>2</sub> ·6H <sub>2</sub> O	25 mg	12.5 mg
Store at 4 °C		

Updated on 11/10/2021

#### AA Amino Acids 100X solution

	For 1 L	For 500 ml
L-glutamine	87.6 g	43.8 g
L-aspartic acid	26.6 g	13.3 g
L-arginine	17.4 g	8.7 g
Glycine	750 mg	375 mg
Store at 4 °C; heat to dissolve		

#### MS Vitamins 100X solution

	For 1 L	For 500 ml
Myo-inositol	10 g	5 g
Nicotinic acid	50 mg	25 mg
Pyridoxine-HCl	50 mg	25 mg
Thiamine-HCl (10 mg/ml)	1000 µL	500 µL
Glycine (100 mg/ml)	20 µL	10 µL
Store at 4 °C		

#### Fe<sub>2</sub>EDTA (Iron Chelate) 100 X Solution

- Dissolve 2.78 g of FeSO<sub>4</sub>·7H<sub>2</sub>O in 350 ml H<sub>2</sub>O
- Dissolve 3.726 g of disodium EDTA (MW: 372.24) in 350 ml H<sub>2</sub>O. Heating is required
- Combine the solutions and bring up to the final volume of 1 L (clear yellow solution)
- Storage: amber bottle (or covered by foil) to protect it from exposure to light; store at 4 °C



Updated on 11/10/2021

## Additional Links

Acetosyringone's function in plant transformation

<https://bmcbiotechnol.biomedcentral.com/articles/10.1186/s12896-018-0459-5>

pCAMBIA 1300 Vector

[https://www.snapgene.com/resources/plasmid-files/?set=plant\\_vectors&plasmid=pCAMBIA1300](https://www.snapgene.com/resources/plasmid-files/?set=plant_vectors&plasmid=pCAMBIA1300)

T-DNA Binary Vectors

<http://www.plantphysiol.org/content/146/2/325>